
PRODUCT INFORMATION

Product Name: T7 Endonuclease I

Code: DS715 **50 rxn** *Research Use Only*

Kit Contents:

Component	Quantity	Storage
①10× T7 Endonuclease I	1 Tube	-20°C
②10× Digestion Buffer	1 Tube	-20°C
③Digestion Control	1 Tube	-20°C

Storage condition:

Store at -20°C

Stable for 12 months from the date of receipt

※③Digestion Control may form sediment. In such cases, use the supernatant after centrifugation.

Introduction:

T7 Endonuclease I is a mismatch DNA cleavage enzyme that detect the editing efficiency of genomic DNA. Loci amplicons where the gene-specific double-strand breaks occur is denatured and reannealed to generate mismatches. The mismatches are cleaved by T7 Endonuclease I and then the resultant bands are analyzed by gel electrophoresis.

Required materials not included with the kit:

- DNase free water
- Block heater
- Tubes, Micro Pipettes & Tips
- Agarose gel, Electrophoresis system & Analysis software

PRODUCT INFORMATION

Methods

Cleavage assay:

1. Combine following samples for re-annealing.


	Sample	Digestion Control
PCR Product	1–8 μ L	1 μ L
10 $\times$ Digestion Buffer	1 μ L	1 μ L
Water	0–7 μ L	7 μ L
Total	9 μ L	9 μ L

2. Place the PCR tube from above in a thermal cycler with a heated-lid and run the following program.

Stage	Temp	Time
1	95°C	5 min
2	95–85°C	–2°C/sec
3	85–25°C	–0.1°C/sec
4	4°C	Hold

3. Add 1 μ L 10 $\times$ T7 Endonuclease I to test samples and water to the control, mix well.
4. Incubate at 37°C for 15 min. Place samples at 4°C immediately to reduce the overreaction.
5. Run agarose gel electrophoresis.
(If not analyzed immediately, add 1 μ L of 100 μ M EDTA to each reaction and store samples at –20°C.)

Cleavage of control:



3% Agarose Gel Electrophoresis